




ProfilerPro Glycan Profiling Quick Guide

Sample Preparation

Denature

1. Thaw and spin Denaturing Plate at 1200g for 1 minute.
2. Carefully remove plate seal.
3. Add 8 μ L of sample (monoclonal antibody) with concentration range of 1.25mg/mL to 7.5mg/mL (10 μ g to 60 μ g total protein) to Denaturing Plate. Mix by pipetting up and down or with a plate shaker.
4. Seal plate carefully with an adhesive plate seal.
5. Spin plate at 1200g for 1 minute.
6. Incubate for 10 minutes at 70C using a PCR machine or heat block.



Digestion

1. Thaw PNGase F Plate.
2. Spin both PNGase F Plate and Denaturing Plate at 1200g for 1 minute.
3. Carefully remove plate seals.
4. Transfer all denatured sample to PNGase F Plate. Mix by pipetting up and down or with a plate shaker.
5. Seal plate carefully with an adhesive plate seal.
6. Spin PNGase F Plate at 1200g for 1 minute.
7. Incubate for 1 hour at 37C using a PCR machine or heat block.

Labeling

1. Thaw Dye Plate.
2. Spin both Dye Plate and PNGase F Plate at 1200g for 1 minute.
3. Carefully remove plate seals.
4. Transfer 8 μ L of Digested sample to the Dye Plate. Mix by pipetting up and down or with a plate shaker.
5. Spin Dye Plate at 1200g for 1 minute.
6. Incubate the unsealed plate for 2 hours at 55C, or until dry using a PCR machine (lid open) or heat block.

Reconstitution

1. Add 100 μ L of molecular grade water to dried samples.
2. Seal plate carefully with an adhesive plate seal.
3. Mix samples on a plate shaker at maximum speed for at least 1 minute.
4. Spin plate at 1200g for 1 minute.
5. Carefully remove plate seal.
6. Run plate on LabChip GXII.

For full Profiler Pro Glycan Profiling User Guide, go to www.caliperls.com/chip-datasheets

ProfilerPro Glycan Profiling Quick Guide

Chip Preparation Procedures

Note: Allow the chip and Gel Matrix ● to equilibrate to room temperature for about 20-30 minutes before use.

Preparing the Buffer Tube

1. Add 750 μ L of molecular grade water to the 0.75mL Buffer Tube.
2. Insert the Buffer Tube into the buffer slot on the LabChip GXII instrument.



Buffer Tube

Preparing the Ladder Tube


1. Add 145 μ L of Ladder Diluent ● to one of the Ladder ● tubes.
2. Vortex at highest speed for about 30 seconds and spin down.
3. Transfer 120 μ L of prepared ladder to the 0.2mL Ladder Tube.
4. Insert the Ladder Tube into the ladder slot on the LabChip GXII instrument.



Ladder Tube

Preparing the Chip

1. Remove reagents from all wells of the chip using a vacuum.
2. Rinse and aspirate all wells with molecular grade water. Repeat, and be careful to aspirate all the water in the wells and any water that may have spilled onto the outside of the chip.
3. Add 75 μ L of Gel Matrix ● to chip wells 2, 3, 7, 8 and 9 and 120 μ L in well 10 using a Reverse Pipetting Technique.
4. Prepare marker solution by adding 125 μ L of Marker Diluent ○ to one of the Marker ● tubes. Vortex at highest speed for about 30 seconds and spin down. Transfer 100 μ L of prepared marker solution to chip well 4. (*Note: Prepare marker solution just before loading the chip in the LabChip GXII and starting the assay. Do not prepare marker solution in advance as the marker signal degrades over time.*)
5. Place the chip in the LabChip GXII instrument to begin the assay.



For full ProfilerPro Glycan Profiling User Guide, go to www.caliperls.com/chip-datasheets