



Accessories for the **RapidTrace** | SPE Workstation

Caliper Life Sciences is proud to offer accessories for the RapidTrace SPE Workstation. For more technical information, pricing and ordering contact your local Caliper Life Sciences Representative by calling 508-435-9500, Fax: 508-435-3439 or email us at: solutions@CaliperLS.com

Racks

There are three racks available. Each RapidTrace requires one of these racks to run. Please select from one of the rack options below.



	DESCRIPTION	QUANTITY	PART NUMBER
Standard Rack	13x100 mm sample tube	1 Rack	50974
	12x75 mm fraction tube	5 Racks	50976
Optional Rack	16x100 mm sample tube 16x100 mm fraction tube	1 Rack	56786
Chilled Rack	13x100 mm sample tube 12x75 mm fraction tube (9 Position Rack)	1 Rack	56536
Magnet Kit	Replacement kit of 10 magnets for encoding racks	1 Kit	54405

Vial Adapters

The vial adapters are designed for use with the standard rack (13x100 mm sample tube, 12x75 mm fraction tube), to collect the fraction directly into a vial rather than a test tube.

Transferring the sample from a test tube to a vial takes additional time and can result in a loss of some of the sample. There are two adapters available, one for an 11 mm vial and one for the 4mL vials.



DESCRIPTION	QUANTITY	PART NUMBER
11mm Vial Adapter	10 Adapters	55091
4mL Vial Adapter	10 Adapters	55827

Empty SPE Columns and Spacers

Empty SPE columns and spacers are used to perform the RapidTrace Cleanup procedure. Each RapidTrace Module comes with ten 3 mL SPE columns and spacers, and two 1 mL SPE columns, and spacers. Additional columns and spacers are available.

DESCRIPTION	QUANTITY	PART NUMBER
Kit of 10 empty 3mL SPE syringe barrel Columns and Spacers	10 Columns and Spacers	59812
Kit of 10 empty 1mL SPE syringe barrel Columns and Spacers	10 Adapters	57145



Test Tubes

Disposable borosilicate glass test tubes are available from Caliper Life Sciences for the RapidTrace.

DESCRIPTION	QUANTITY	PART NUMBER
12 x 75mm Test Tubes, uncapped	1,000/case	44651
13 x 100mm Test Tubes, uncapped	1,000/case	40707
16 x 100mm Test Tubes, uncapped	1,000/case	40708



O-Rings

The O-Rings used on the rapidTrace are subject to wear and will need to be replaced occasionally.



DESCRIPTION	QUANTITY	PART NUMBER
3mL SPE O-Ring, CHEMRAZ™	1 O-Ring	45096
1mL SPE O-Ring, CHEMRAZ	1 O-Ring	60969
Cannula O-Ring, CHEMRAZ	1 O-Ring	45096
Waste Station O-Ring, CHEMRAZ (3 in a Waste Station)	1 O-Ring	58526

Documentation

DESCRIPTION	PART NUMBER
RapidTrace SPE Workstation Installation and Quick Reference Manual	55078
RapidTrace Validation Documentation Set	56846

Miscellaneous

DESCRIPTION	QUANTITY	PART NUMBER
3mL SPE Gel Adapter Kit	Quantity 1	58728
1mL SPE Gel Adapter Kit	Quantity 1	58759
Gas Conversion Lines	1 Kit	58262
1mL Adapter Sleeves	Kit of 11	59807
Tubing Weight Kit for Reagent Lines	Set of 8 weights	50973
Check Valve Kit for Manifold Reagent Lines	Set of 8 Check Valves	55076
Software Rev 1.20 Upgrade Kit	Quantity 1	58777



Corporate Headquarters

68 Elm Street
Hopkinton, MA 01748-1668

Tel: 1.508.435.9500 / 1.877.LabChip

Fax: 1.508-435-3439

Email: cust.support@caliperLS.com

www.caliperLS.com

©2008 Caliper Life Sciences. All rights reserved.
Caliper, RapidTrace and TurboVap are tradenames
and/or trademarks of Caliper Life Sciences, Inc.

RT-BR-02 Apr08