

# BioRobotix™ 96 Pipet Tips

The Official Tip For Use With Caliper Life Sciences Instruments

## Caliper Automation Certified Pipet Tips

Caliper Life Sciences is proud to be partnered with Molecular BioProducts to provide Caliper Automation Certified Pipet Tips under the BioRobotix brand. 'Caliper Automation Certified' is the exclusive trademark of Caliper Life Sciences. This seal means that you are getting the very best product available.



- Guaranteed box-to-box uniformity and trouble-free performance
- Nested tray design is three times more space efficient compared to traditional racked tips!
- Microplate shaped trays facilitate robotic handling.
- Available in racked or nested formats
- DNase and RNase free.

## Caliper Life Sciences has introduced a new & improved tip offering for the Sciclone ALH HV. These new tips bring expanded functionality including:

- Nested 30-Stacks to fit your Twister II Storage Design.
- SBS standardized footprint to allow use of tip positions interchangeably with microplates.
- 50ul volume tips for applications requiring small volume liquid handling precision.
- Wide bore sterilized tips for cell based assay or DNA handling.

Why not schedule tip deliveries to your lab? Our supply contracts offer storage efficiency, convenience, and additional savings! Ask about our supply contract programs.



# BioRobotix 96 Pipet Tips

## Specifications

Part Number	Description	Units Format
109081	200ul 96 Rack Tips	10 Racks of 96 Tips
109082	200ul Sterile 96 Rack Tips	10 Racks of 96 Tips
111623	200ul Wide Bore, Sterile 96 Rack Tips	10 Racks of 96 Tips
111426	150ul Barrier Sterile 96 Rack Tips	10 Racks of 96 Tips
109079	100ul 96 Rack Tips	10 Racks of 96 Tips
109080	100ul Sterile 96 Rack Tips	10 Racks of 96 Tips
111624	80ul Barrier Sterile 96 Rack Tips	10 Racks of 96 Tips
111625	100ul Carbon 96 Rack Tips	10 Racks of 96 Tips
105641	50ul 96 Rack Tips	10 Racks of 96 Tips
105642	50ul Sterile 96 Rack Tips	10 Racks of 96 Tips
105645	35ul Barrier Sterile 96 Rack Tips	10 Racks of 96 Tips

Part Number	Product	Units per Sleeve
105649	200ul Nested Twister Tray 96 Tips	30 Trays of 96 Tips
105650	200ul Sterilized Nested Twister Tray 96 Tips	30 Trays of 96 Tips
105647	100ul Nested Twister Tray 96 Tips	30 Trays of 96 Tips
105648	100ul Sterilized Nested Twister Tray 96 Tips	30 Trays of 96 Tips
105643	50ul Nested Twister Tray 96 Tips	30 Trays of 96 Tips
105644	50ul Sterilized Nested Twister Tray 96 Tips	30 Trays of 96 Tips

