



Tempo™ 384 Pipet Tips

The Official Tip For Use With Caliper Life Sciences Instruments

Caliper Automation Certified Pipet Tips

are now available in 25ul volumes.

'Caliper Automation Certified' is the exclusive trademark of Caliper Life Sciences. This seal means that you are getting the very best product available.



- Guaranteed box-to-box uniformity and trouble-free performance
- Nested tray design is three times more space efficient compared to traditional racked tips!
- Microplate shaped trays facilitate robotic handling.
- Available in racked or nested formats
- DNase and RNase free.

Caliper Life Sciences has introduced a new & improved 25ul tip offering for the Sciclone ALH LV. These new tips bring a number of added benefits including:

- Attach rows, columns, or individual tips -- The LV head can now be used to carry out both row-wise and column-wise serial dilutions.
- Extract 96 tips from 384 rack or tray - Conserves storage space, and minimizes product SKU requirements.
- Access 1536 well plates - Confidently perform 384 to 1536 transfers.
- 0.5ul Transfer Volumes - Volume transfer ranges are the same for tips or cannulae.

Caliper Automation Certified Tempo 384 Pipet Tips Include:

Part Number	Description	Unit Format
108849	25ul Pipet Tips, Racked	10 Racks of 384 Tips
108851	25ul Pipet Tips, Sterilized, Racked	10 Racks of 384 Tips
108848	25ul Pipet Tips, Nested	35 Trays of 384 Tips
108850	25ul Pipet Tips, Sterilized, Nested	35 Trays of 384 Tips

Why not schedule tip deliveries to your lab? Our supply contracts offer storage efficiency, convenience, and additional savings! Ask about our supply contract programs.

Caliper Life Sciences, Inc. - 68 Elm Street - Hopkinton, MA 01748-1668 - TEL: 508-435-9500 - FAX: 508-435-3439 - Website: www.CaliperLS.com

©2004 Caliper Life Sciences, Inc. Caliper, Automation Certified and LabChip are registered trademarks, and the Caliper Logo and the phrase "Working Innovation" are trademarks of Caliper Life Sciences, Inc..

Con-DS-01 06/04

